

# Roadster Deluxe



**Weight Capacity:** 260 lbs

**Max. Speed:** 4mph

**Range:** 6mi.

**Turning Radius:** 35.4"

**Ground Clearance:** 2"

**Controller:** PG S-Drive 45

**Battery:** 15Ah

**Charger:** 1.8A off-board

**Gradient:** 6°

**Front Wheel:** 8"(200x50) PU tire

**Rear Wheel:** 8"(200x50) PU tire

**Motor:** DC 24V, 110W(Rating)/420W(Max)

**Brake:** intelligent, regenerative, electromagnetic brakes



seat dimensions inches			seat height inches		wheel base inches		
D	E	F	G	H	I		
18	15	15	19-21	16 -18	31		
overall dimensions				weight / pounds			
A	B	C	base	battery/ea	seat	ttl w/batt.	capacity
40	21	35	64	11	15	99	260